

Scio Four Family Gas Measurement Modules

O₂, CO₂, N₂O and volatile anaesthetic agents at a glance: Scio Four Family can be used with an Infinity monitor anywhere you need it.



Dräger. Technology for Life®

Benefits

Accurate real-time gas measurement

Whether inspiratory or expiratory values for O₂, CO₂, N₂O and anaesthetic agents: Dräger Scio Four Family delivers precise and reliable values which are transferred to a Dräger Infinity monitor and displayed in real time.

Scalable solutions for every application

The Scio Four Family consists of four different versions: options include automatic anaesthetic agent identification as well as agent mixture measurement and integrated O_2 measurement (see technical data on page 5).

The modular design of the Scio Four Family makes it both flexible and scalable. If you have an integrated O_2 measurement in your Fabius device, a gas measurement module without oxygen sensor will meet your needs. This means you get just the system you need – no more, no less.

Enhanced intraoperative patient safety

The Scio Four Family does more than just measure gas concentrations in the breathing circuit. In combination with Dräger Infinity monitors, MAC values for anaesthetic agents are displayed, helping you to avoid underand overdosage. This means you can maintain optimal gas concentrations while making sure that ventilation is adequate.

You also have the option of measuring inspiratory and expiratory O_2 concentrations using a paramagnetic oxygen sensor. This consumption-free measurement method – unlike electrochemical sensors – results in reduced costs for you.

Waterlock® 2: enhanced contamination protection

The Scio Four Family comes equipped with the WaterLock 2 water trap. The 2 µm pores of its dual hydrophobic membrane design let gases pass freely while effectively preventing microbial contamination and condensation. Even under challenging conditions such as minimal flow anaesthesia, where heavy condensation can be a problem, accurate capnography is assured. Sample and purge gases remain free from condensation. Highly economical, the WaterLock 2 can be emptied at any time and only requires monthly replacement.

Made for Infinity® and beyond

The Scio Four Family is specially designed for use with Dräger Infinity monitoring systems^{*}. Just plug in the gas measurement module and your Infinity monitor will display all the values generated. Your Dräger representative can help you determine which Dräger Infinity monitors are best for you.

* e.g. Delta, Delta XL and Vista XL (not available in all countries)

System components



Dräger Fabius® GS premium

The Dräger Fabius[®] GS premium is an anaesthesia workstation that is simple to use, highly efficient and ready for the future. It features a solid design with modular architecture plus a wide range of ventilation capabilities. Your Fabius[®] GS premium can be customised to suit your exact needs.

Dräger Fabius® Tiro

Get the most out of even the smallest spaces with a compact yet fully featured anaesthesia solution designed for use in a variety of specialised environments.



D-6838-2014

Infinity[®] Delta XL

With a 12.2" (310 mm) colour screen, the Delta XL monitor can continuously monitor patients both at the bedside and during transport – thus eliminating the need for separate transport monitors. Supports all patient types and acuity levels throughout the hospital.

Didger	Infinity Vista
petriladadadad	ruhruhruh 📑 🛛 😽 📒 🖯
fitriladadadad	rdride 🖉 🗄
mm	MM 14
	15
38 20*	15*
00	C

Infinity[®] Vista XL

Offering true Dräger quality, the Vista XL provides the capabilities you need to monitor adult, paediatric and neonatal patients. It's easy to learn, easy to use and can function as a stand-alone monitor or a networked device to maximise your connection to vital patient information.

Accessories



WaterLock[®] 2

Perfect protection for precise gas measurement. Dräger WaterLock[®] 2 safely stops water from getting into the multi-gas sensor. The membrane technology developed by Dräger for the WaterLock[®] 2 stops any bacteria or germs from getting into the gas measurement system. The WaterLock[®] 2 is also safe and simple to empty – with a further advantage in handling and hygiene.

Related Products



Vamos Family

The expandable Vamos Family lets you choose your gas monitoring solution today while keeping your options open for the future. In addition to CO_2 and SpO_2 your monitoring capabilities can be expanded to include all five common anaesthetic agents and N_2O .

Technical Data

Operational characteristics

Weight		2.9 kg (6.4 lbs) w/o O₂ sensor		
		3.2 kg (7.0 lbs) w O ₂ sensor		
Dimensions (H x W x D)		11.5 cm x 19.0 cm x 27.0 cm (4.8 in x 7.5 in x 10.5 in)		
Power supply		100 to 240 V ~ 50/60 Hz		
Device noise emission		≤45 dB(A)		
Device interfaces		1x RS 232 (MEDIBUS protocol)		
Ambient conditions				
Temperature		10 to 40°C (50) to 104°F)	
Ambient pressure		700 to 1,100 hPa (10.2 to 15.9 psi)		
		Equivalent to 3,000 m (9,843 ft) height		
Measurement specification				
O ₂	Range		0 to 100 vol.%	
	Accuracy		±(2.5 vol.% + 2.5% rel.)	
	Rise time (t _{10/90})		<500 ms	
	Range		0 to 10 vol.%	
	Accuracy		±(0.43 vol.% + 8% rel.)	
	Rise time (t _{10/90})		<300 ms	
N ₂ O	Range		0 to 100 vol.%	
	Accuracy		±(2 vol.% + 8% rel.)	
	Rise time (t _{10/90})		<300 ms	
Anaesthetic agents	Range		0 to 8.5 vol.% (isoflurane and halothane)	
			0 to 10 vol.% (sevoflurane and enflurane)	
			0 to 20 vol.% (desflurane)	
	Accuracy		±(0.20 vol.% + 15% rel.)	
	Rise time (t _{10/90})		<350 ms	
			<450 ms (halothane)	
Respiratory rate	Range		0 to 100/min	
	Accuracy		0 to 60/min: ±1/min	
			>60/min: not specified	
Device classification				
Protection class (IEC 60601-1)		Class 1, type E	continuous operation	
Ingress of liquids and particulate	ress of liquids and particulate matter (IEC60601-1)		IP21	
Classification in accordance with EU		lla		
Directive 93/42/EEC Appendix IX	X			
UMDNS code / GMDN code		17-445 / 41688	3	
Device name	Automatic agent II	D	Integrated O ₂ measurement	
Scio Four Oxi plus	yes		yes	
Scio Four plus	yes		no	
	yes			
Scio Four Oxi	no		yes	

Notes

compatible | Subject to modifications | © 2015 Drägerwerk AG & Co. KGaA tally free iications & Sales Marketing | PP | LE | Printed in Germany | Chlorine 90 66 495 | 15-09-3 | Comr

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA Moislinger Allee 53–55 23558 Lübeck, Germany www.draeger.com

Manufacturer:

Drägerwerk AG & Co. KGaA Moislinger Allee 53-55 23558 Lübeck, Germany

Locate your Regional Sales Representative at: www.draeger.com/contact



REGION EUROPE CENTRAL AND EUROPE NORTH Drägerwerk AG & Co. KGaA Moislinger Allee 53-55 23558 Lübeck, Germany Tel +49 451 882 0 Fax +49 451 882 2080 info@draeger.com

REGION EUROPE SOUTH Dräger Médical S.A.S. Parc de Haute Technologie d'Antony 2 25, rue Georges Besse 92182 Antony Cedex, France Tel +33 1 46 11 56 00 Fax +33 1 40 96 97 20 dlmfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA Branch Office P.O. Box 505108 Dubai, United Arab Emirates Tel +971 4 4294 600 Fax +971 4 4294 699 contactuae@draeger.com

REGION ASIA / PACIFIC

Draeger Medical South East Asia Pte Ltd. 25 International Business Park #04-27/29 German Centre Singapore 609916, Singapore Tel +65 6572 4388 Fax +65 6572 4399 asia.pacific@draeger.com

REGION CENTRAL

AND SOUTH AMERICA Dräger Panama Comercial S. de R.L. Complejo Business Park, V tower, 10th floor Panama City Tel +507 377 9100 Fax +507 377 9130 contactcsa@draeger.com

CANADA

Draeger Medical Canada Inc. 2425 Skymark Avenue, Unit 1 Mississauga, Ontario, L4W 4Y6 Tel +1 905 212 6600 Toll-free +1 866 343 2273 Fax +1 905 212 6601 Canada.support@draeger.com