



Reliable, every day

The Primus 1058 Life was designed with a major focus on functional reliability and economic efficiency - for today and tomorrow. Demonstrate your ability to identify excellent value for your money and join thousands of dentists throughout the world who report that their Primus 1058 "runs and runs and keeps on running".

Highlights at a glance:

The patient chair

1. More streamlined chair design in modern dental white
2. Maximum load up to 185 kg
3. Integrated Trendelenburg movement for patient comfort
4. Extended highest/lowest positions
5. Optimised comfort from the soft cushion shape, optional RELAXline soft upholstery

The dentist and assistant elements

6. Dentist element in a modern design
7. Simple operating concept with more direct selection buttons
8. Assistant element includes clearly structured film keypad and improved operating comfort

Patient communication

9. KaVo ProXam iCam intraoral camera or KaVo DIAGNOcam Vision Full HD*
10. New high-resolution KaVo screens, 19" and 22"
11. The new CONEXIO system software combines patient management software, a 2D radiograph archive and live intraoral images in a unique operating concept.

*Products not shown

Equipment options

	TM	S	Cart		TM	S	Cart
Patient chair							
Standard	●	●	●		●	●	●
COMPACTchair	○	○	○		●	●	●
Soft upholstery	●	●	●		○	○	○
RELAXline Soft upholstery	○	○	○		○	○	○
Backrest							
Progress	●	●	●		○	○	○
Comfort ²	○	○	○				
Headrest							
Headrest	○	○	○				
Double-jointed with rotary wheel	●	●	●				
Double-jointed with push button	○	○	○				
Armrest							
Right	○	○	○				
Left	○	○	○				
Dentist element							
Turbine hose	○	○	○				
KL-Motor 701	○	○	○				
KL-Motor 703 LED™	○	○	○				
Triple function handpiece	○	○	○				
Minilight triple function handpiece	○	○	○				
Multifunctional handpiece	○	○	○				
COMFORTdrive™	○	○	○				
SMARTdrive™	●	●	●				
PiezoLED™ ultrasonic scaler	○	○	○				
SONICflex™ sonic scaler	○	○	○				
MemoSPEED display indication-rel.	○	○	○				
Panoramic X-ray viewer ⁴	○	○	○				
(mounted to light mounting post)							
Panoramic X-ray viewer	-	-	○ ³				
Tray holder for US tray	x	x	○				
Tray holder for single standard tray	x	-	x				
Tray holder for two standard trays	○	○	-				
Spray heating for instruments	○	○	○				
Bell function key	○	○	○				
Patient communication							
KaVo Screen One	○	○	○				
KaVo Screen HD	○	○	○				
CONEXIO™	○	○	○				
KaVo ProXam iCam	○	○	○				
KaVo DIAGNOcam Vision Full HD	○	○	○				
Assistant element							
Spray mist suction	●	●	●		●	●	●
Saliva ejector	●	●	●		○	○	○
2nd saliva ejector	○	○	○		○	○	○
Triple function handpiece	○	○	○		○	○	○
Multifunctional handpiece	○	○	○		○	○	○
Minilight triple function handpiece	○	○	○		○	○	○
Curing light Satelec MiniLED	○	○	○		○	○	○
Support for tray holder	○	○	○		●	●	●
Right/left swivelling and height-adjustable assistant element	●	●	●		○	○	○
Right/left swivelling and height-adjustable assistant element	○	○	○				
Suction tube guide	○	○	○				
Selective support kit	○	○	○				
Unit body							
Porcelain spittoon bowl	x	x	x				
Glass spittoon bowl	x	x	x				
VACUstopp	●	●	●				
Permanent germ reduction ⁵	●	●	●		○	○	○
Intensive germ reduction ⁶	○	○	○		●	●	●
HYDROclean	○	○	○				
Warm water heater	○	○	○				
Supply system							
Water block DVGW	x	x	x				
Water block compact	x	x	x				
Water bottle DVGW	x	x	x				
Low-pressure regulator	○	○	○				
Connection to external device	○	○	○				
Disposal system							
External suction	x	x	x				
DÜRR amalgam separator	x	x	x				
DÜRR automatic separator	x	x	x				
Solids collector kit	x	x	x				
Miscellaneous							
KaVo PHYSIO Evo/ Evo F	○	○	○				
KaVo PHYSIO One	○	○	○				
Operating light KaVoLUX™ 540 LED	○	○	○				
Operating light MAIA LED	○	○	○				
Light attachment bar	○	○	○				
Adapter for light attachment bar	○	○	○				
Patient service box	○	○	○				
Service table 1568	-	-	○				
DEKASEPTOL Gel basic set	●	●	●				

● Standard equipment ○ Optional equipment x Must be selected - Not available

1. Not with COMPACTchair patient chair
2. Only with COMPACTchair patient chair
3. Not suitable for placing under a tabletop

4. Not in conjunction with multimedia
5. Only with DVGW water block or DVGW water bottle
6. Only in conjunction with DVGW water block